

High Pressure Oxygen Concentrator(0.4Mpa)



Specifications

Flow Rate	0-20L/min
Purity	93±3%
Weight	148kg
Working Voltage	230V
Power Consumption	2200W
Outlet Pressure	0.4Mpa
Noise Level	62dB
Dimension	800*700*1030mm
Dehumidifier	Cooler-dehumidifier
Emergency Signals and Alarms	Power failure alarm, Low purity alarm, Low pressure alarm, Over Heat alarm, High pressure malfunction alarm