

# DMH Series



9kg light weight, portable and convenient



Oxygen flow adjustable from 1L-5L/min



Digital flowmeter adjustable by 0.1L/min interval



Super low noise for a quiet experience



Advanced oxygen purity real-time display



## Specifications

Flow Rate	0.5-5L/min , rated flow:2L/min
Maximum	5L/min
Net Weight	9KG
Purity	93±3%
Power Consumption	242W
Outlet Pressure	0.04±10%
Noise Level	≤42dB
Dimension	344*230*300mm
Color	Accept customization
OEM&ODM	Acceptable
Standard Function	Over heat alarm. Power failure alarm.Timing function.Working hours display,Low purity alarm , Purity real-time display
Optional Function	Nebulizer Function