

# VH5-B Series



14kg total weight -easy lift and convenient move



Digital flowmeter adjustable by 0.1L/min interval



SpO2 sensor or nebulizer function available for choosing



Reliable user experience from low purity alarm



## Specifications

Flow Rate	0.5-5L/min
Net Weight	14.3kg
Purity	93±3%
Working Voltage	230V
Power Consumption	285W
Outlet Pressure	0.04-0.08Mpa
Noise Level	≤48dB
Dimension	360*260*585mm
Color	Blue (Accept customization)
OEM&ODM	Acceptable
Standard Function	Low purity alarm, Pressure failure alarm,Over heat alarm,Power failure alarm,Timing function
Optional Function	Nebulizer function,SpO2