

VH5 Series



Digital flowmeter adjustable by 0.1L/min interval



Large LCD screen display



Advanced oxygen purity real-time display



SpO2 sensor and nebulizer function available for choosing



Specifications

Flow Rate	0.5-6L/min , rated flow:5L/min
Purity	93±3%
Working Voltage	230V
Power Consumption	480W
Net Weight	18KG
Outlet Pressure	0.04-0.08Mpa
Noise Level	≤48dB
Dimension	375*252*574mm
Color	Grey (Accept customization)
OEM&ODM	Acceptable
Standard Function	Low purity alarm, Over heat alarm, Power failure alarm, Timing function, Working hours display, Pressure failure alarm,Purity real-time display
Optional Function	Nebulizer Function , SpO2