



OXYGEN CONCENTRATOR MEDICAL

V Series

WHO and UNICEF selected models with high reputation



Intelligent self-checking system with error code displayed on LCD



Super low noise for a quiet experience



Option for 24 hours/365 days continuously working



Sustainable purity and flow rate for the device lifespan



Specifications

Flow Rate	0-5L/min	0-8L/min	0-10L/min
Net Weight	21/24kg	24kg	25kg
Purity	93±3%	93±3%	93±3%
Working Voltage	230V	230V	230V
Power Consumption	350W/480W	650W	650W
Outlet Pressure	0.04-0.08Mpa	0.04-0.08Mpa	0.04-0.08Mpa
Noise Level	45dB/≤50dB	≤50dB	≤50dB
Dimension	390*337*620mm	390*337*620mm	390*337*620mm
Color	Dark Grey (Accept customization)		
OEM&ODM	Acceptable		
Standard Function	Over heat alarm, Power failure alarm, Pressure failure alarm, Timing function, Working hours display		
Optional Function	Low purity alarm, Nebulizer function, Dual flow		