



HYPERBARIC CHAMBER FOR SPORT



Hard Type Hyperbaric Chamber



Product Features

- Three-in-one control unit that combines oxygen concentrator, air compressor, and dehumidifier together has one-click operation function
- Plasma disinfection (high efficiency and low consumption)
- Provide internal/external large LCD screen, allow you to operate the equipment by yourself, Easy to control the device inside the chamber
- Multiple safety control
- Timing Function (make sure the security of oxygen session)
- Automatic decompression (balance pressure in chamber)
- Manual decompression for acceleration. (for emergency)
- Regulate oxygen flow (regulate oxygen flow at any time in your oxygen session)
- LCD Touchscreen, easy to operate

Specification --- Control unit

Type	Oxygen Flow	Air Flow	Size(cm)	Weight	Voltage	Power Consumption	Display screen
3-in-1 control (oxygen concentrator, compressor, dehumidifier)	10L/min	140L/min(1.3ATA) 132L/min(1.5ATA)	L57*W46*H70	78kg	100~240V	1100W	LCD touchable

Specification --- Chamber

Type	Pressure	Weight	Size(CM)	Material	Color	Pressure Gauge	Frame
Soft chamber (lying type)	1.3ATA 1.5ATA	130kg 150kg	Ø75*L216	Aluminium+ polycarbonate	White	Inside and outside	LCD touchable, lighting, wireless communication system, Plasma disinfection, Safety pressure relief, One button start