

# Oxygen Filling System



This oxygen filling system is available for medical, homecare, industrial, and other scenarios where delivery of oxygen cylinders is difficult. It is suitable for EMS oxygen supply for life rescue in disaster areas as well as for health care in low purity underground or plateau.

## Features

1. Dry, clean, and oil-free
2. Energy-efficient and Cost-effective
3. Small size and light weight, easy to operate
4. Efficient filling: 3 hours to fill a 10-liter oxygen cylinder
5. High compression ratio up to 138bar (2000psi/13.8Mpa)
6. Separately use or expand into a larger system to meet the demand for a large amount of oxygen consumption