

# Multiplace Hard Type Hyperbaric Chamber



## Product Features

- Three-in-one control unit that combines oxygen concentrator, air compressor, and cooler dryer, has one-click operation function
- Plasma disinfection (high efficiency and low consumption)
- Provide internal/external large LCD screen, allow you to operate the equipment by yourself. Easy to control the device inside the chamber
- Multiple safety control
- Timing Function (Make sure the security of oxygen session)
- Automatic decompression (balance pressure in chamber)
- Manual decompression for acceleration. (For emergency)
- Regulate oxygen flow (Regulate oxygen flow at any time in your oxygen session)
- LCD screen, easy to operate
- Oxygen concentration detection in chamber

## Specification -- Control unit

Type	Oxygen Flow	Air Flow	Size(cm)	Weight	Voltage	Power Consumption	Display screen
3-in-1 control (oxygen concentrator, compressor dehumidifier)	10L/min	85L/min	L90*W80*H105	200kg	100~240V	2200W	LCD touchable

## Specification -- Chamber

Type	Pressure	Weight	Size(CM)	Material	Color	Pressure Gauge	Frame
Hard chamber	1.3ATA	880kg	L200*W124*H200	Carbon steel+ polycarbonate	pearl grey + Light grey	inside and outside	LCD touchable, Lighting, Interphone communication system, Plasma disinfection, Safety pressure relief, One button start