

# Portable Oxygen Concentrator

-  light weight design of only 1.95kg
-  Nearly 4 hours long duration of battery
-  Both pulse and continuous working mode
-  1-5 gears flow regulation



## Specifications

Noise Level	55dB(A)
Purity	≥90%
Continuous Flow	1L/min
Pulse Flow Rate	Pulse Flow Setting 1-5
Max Pressure of Outlet Port	40Kpa±10%
AC Power	AC input: 100-240V~50-60Hz
DC Power	DC input: DC19V 6.31A 120W
Rechargeable Battery	Voltage: 14.8V 6.4Ah 95Wh
Battery Duration	Up to 4 hours
Net Weight	1.95kg
Dimension	217*85*156mm

**PORTABLE OXYGEN  
CONCENTRATOR**